



# Clover 60 Vet

Affordable Excellence for All

**wisonic**



# Vet Application

Pet Hospitals and Clinics

1. General application
2. Sterilization

Horse Riding Club

1. Physical examination
2. Breeding and monitoring equine pregnancy
3. Tendon scanning

Livestock

1. Bovine, equine, ovine pregnancy
2. Musculoskeletal
3. Meat and fat quality rating

Aquaculture

1. Gender screening
2. Breeding selection

Safari & Zoo

1. Regular physical examination
2. Reproduction



# Dedicated Transducers



Convex C5-1  
Medium to large animals  
Abdomen, reproductive system



Linear L15-4  
MSK of medium to large animals  
Small animals whole body



Micro-Convex C8-3  
Small animals abdomen and heart



Phased P7-3  
Small to medium animals heart




Linear L9-4V  
Intracavity, reproductive system

# Excellent Image Quality




Cat Stomach



Dog Pancreas



Dog Heart



Cat Heart



Cat Kidney



Cat Small Intestine

# Broad-Range of Transducers

Clover is capable of performing versatile applications from small animals to large animals with excellent image quality.

